



## Station Security Door



> Pevco pneumatic tube **station security doors** offer added protection for high-value blood products and pharmaceuticals.

- **Restrict unauthorized access to pharmaceuticals**
- **Ensure patient privacy at vulnerable public area stations**
- **Protect against unauthorized station use**

Pevco security doors can be added to existing stations or installed on new stations.

When additional security is needed on a Swisslog™ or other non-Pevco station, a stand-alone door is available which operates independently of existing station controls through a door-mounted 10-digit keypad.

*Swisslog™ pneumatic tube stations are manufactured by Swisslog™ which has no affiliation with Pevco.™*

With a Pevco security door in place, hospital staff must enter an authorization code to release the durable electromagnetic lock to access the station's contents. Once unlocked, the door stays open, leaving hands free to remove the carrier that has arrived. If the door is not manually closed within 60 seconds, a programmable alarm rings as a reminder to close the door. Shortcut keys can be programmed to speed access or users can open the door through the user options function on the main station control pad.

Made of the same durable, powder-coated steel as Pevco's stations, the security door is built to withstand the most demanding environment. The clear, shatter-proof polycarbonate window is secure yet allows users to easily see when a carrier has arrived. For added security, an alarm is triggered if the door lock is ever forced open.

For more information on **station security doors** and other reliable Pevco products, visit [www.pevco.com](http://www.pevco.com) or call 1.800.296.7382.